

ACTIVITY

The electrostatic properties of this uniquely formulated product attract debris and microorganisms away from the skin like a magnet. The elimination of these stress factors allows the normal and healthy function of skin to be restored.

DIRECTIONS

APPLICATION TO WET COAT: Wash animal at least once a week and more if necessary. Thoroughly wet the affected area with water. Spray SuperDerm generously on the area. Gently massage into the skin and allow to stand for 5 minutes before rinsing well with water. Repeat daily until improvement is observed.

APPLICATION TO DRY COAT: Spray dry coat of affected area generously with SuperDerm until very wet. Let stand and dry. Repeat one to three times daily until improvement is observed.

WOUND CLEAN-UP: Clean and rinse area thoroughly with water. Gently apply SuperDerm with a sponge or cloth. Let stand and dry. Repeat on a daily basis.

IF CONDITION DOES NOT IMPROVE,
CONTACT YOUR VETERINARIAN

WARNING: FOR EXTERNAL USE ONLY



CHOICE of CHAMPIONS
USA LLC®

SUPERDERM SOLUTION®

SKIN RESTORATIVE

Helps restore skin due to stress from dermatitis, hot spots, fungus, mange, abrasions, and other skin intrusions.

FOR ALL CLASSES OF HORSES

Net Contents: Quart (32 fl. oz)

SuperDerm Solution® is formulated through a complex, multi-staged, synthesis process using USDA certified bio-based products to produce a wetting agent with a base of colloidal micelles.

INGREDIENTS:

Processed extracts of coconut, corn, soybean, and tree sap formulated with amino acids, fatty acids, alkanolamines, and nonionic surfactants.

ATTRIBUTES:

- Non-toxic, non-carcinogenic
- Safe for human contact
- Hypoallergenic
- Readily biodegradable
- No dyes or perfumes



STORAGE:

Store in a temperature range from 40-80 degrees Fahrenheit and out of direct sunlight.

SuperDerm Solution® is a trademark of Choice of Champions USA LLC. All rights reserved.

Choice of Champions USA LLC
Wellington, Florida • 800-868-1077
SuperDermSolution.com

